

Intensive Care Patient Monitor

MP1000NT *Plus*

Multi-parameter modules

ECG, SpO₂, NIBP, 2-IBP, 2-TEMP, RESP, EtCO₂

Arrhythmia analysis and ST segment analysis

72 hours tabular and graphical trends


Save up to 30 events management

High Resolution 10.4 inch TFT LCD screen

User configurable page layouts

Connection to MEKICS Central Monitoring System via cable network.





△ Large font display



△ Event manager



△ Dose calculator



△ Control knob and 6 Hot keys



△ Multi-parameter modules



△ Wall mount

SPECIFICATIONS

Features

Display Monitor

10.4" Color TFT LCD screen & Encoder, Up to 7 waves

Functions

ECG, Resp, SpO₂, NIBP, 2IBP, 2Temp, EtCO₂
15 Kinds Arrhythmia, ST, Pace Maker Detect
30 Event with 30 sec wave recording / 72 hrs Trend
Virtual Screen, Color Config, Drug Dose Calc, HL7 support
IBP cable, EtCO₂ module.
Thermal Printer (58 mm thermal 3wave)
Cart for MP1000NTP

Optional

Display Data

Parameters

Setting parameters/Patient status parameters/Alarm status

Graphic waveform

ECG (I, II, III, aVR, aVL, aVF, V), SpO₂, IBP, RESP, EtCO₂

Trend

72 hrs Graphic trend of Vital signs

Event Data

Max. 30-arrhythmia event datas(30 sec recording Wave)

Special Functions

Event Managing

Reviewing Functions for wave data of each Event. Max 30 Events.

Auto Printing

Auto Printing in Alarm status

Drug Dose Calc.

HL7 Support

Electrical

Power Source

AC

100-240 VAC, 50 / 60 Hz

Consumption

Max 150 VA

External DC

Type

+12 V / 4 A DC input (Car Cigar Jack)

Battery

Operation Time

Li-ion Battery (DC 10.8 V, 2200 mAh)
1 hrs (Std)

Communication / Output

RS232

COM1

115200 BPS for CMS

LAN

10 Mbps

Physical

Weight

5 Kg without Printer

Dimensions

W 282 x H 315 x D 165 mm

Corporate Headquarters

MEKICS CO., LTD

21, Sangjiseok-gil

Paju-si, Gyeonggi-do

10911, Korea

T. 82 70 7119 2520

F. 82 70 5052 5800

E. sales@mek-ics.com

www.mek-ics.com

Distributor / 대리점

CE 2460

MEK
Intensive Care System